

MEDIBUS for Dräger Pediatric Devices

Instructions for Use

Contents

For Your Safety and that of Your Patients.....	3
Intended Use.....	4
MEDIBUS Specification for Incubator 8000 IC/SC/NC..	5
Commands.....	6
Measured Data.....	7
Alarm Messages.....	7
Device Settings.....	9
Text Messages.....	9
 MEDIBUS Specification for Babytherm 8004/8010.....	 10
Commands.....	11
Measured Data.....	12
Alarm Messages.....	12
Device Settings.....	14
Text Messages.....	14
 MEDIBUS Specification for Babylog 8000.....	 15
Commands.....	16
Realtime Data.....	17
Measured Data.....	17
Supported Alarm Messages.....	18
Device Settings.....	18
Text Messages.....	19

For Your Safety and that of Your Patients

For correct and effective use of the apparatus and to avoid hazards it is essential to read the following recommendations and to act accordingly:

Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

Liability for proper function or damage

The liability for the proper function of the software protocol is irrevocably transferred to the owner or operator if the software protocol is used in a manner not conforming to its intended use.

Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

Dräger Medizintechnik GmbH

Intended Use

MEDIBUS

MEDIBUS is a software protocol intended to be used by two medical devices for exchanging data and control functions via their RS 232 interfaces.

This part of manual describes device dependent supported commands and data sets, port hardware and configurations for Dräger Pediatric Devices.

For a general description of the protocol please refer to the Instructions for Use "Dräger RS 232 MEDIBUS Protocol Definition" (order-no. 90 28 258).

Any data transmitted via the MEDIBUS interface are intended only for information purposes and should not be used to derive therapeutical decisions.

Data Formats

The underscore character used in the format column in lists of measured data and alarm limits is transmitted as an ASCII "space" character (20H).

A '*' ahead the format indicates that the value may be negative. In that case a '-' character will appear at the first space of the respective format.

Be aware that in a small number of cases the used format for a certain item may be different!

Languages

Alarm phrases and text messages are given in the languages:

D	german	NL	dutch
GB	english	E	spanish
US	US-english	J	japanese
F	french	S	swedish
I	italien		

Alarm Phrases

Some alarm phrases contain abreviations as follows:


ASCII Short Form	Meaning
\$&	LOW
"#	HIGH
'@	ALARM

Alarm Priorities

The alarm source is responsible for the priority.
The priority may be variable even if from the same source. (E.g. the priority may increase the longer an alarm is pending.)

MEDIBUS Specification for Incubator 8000 IC/SC/NC

Device Connection



Port Specification

Connector	
Type	RS-232-C
Pins	9 pin Sub D (female)
	1 Housing
	2 RXD
	3 TXD
	5 GND
Galvanic Isolation	1.5 kV
Location	rear side of Incubator
	8000 IC/SC/NC
	Label: Baby Link®
To connect a PC to Incubator 8000 IC/SC/NC the "Medi-Cable" 83 06 488 is recommended.	

Port Configuration

Baudrate	9600 Baud
Databits	8
Startbits	1
Stopbits	1
Parity	none

Device Identification

ID Number	Name	MEDIBUS-Version
8000	INCUBATOR 8000 SC/NC	03.00
8001	INCUBATOR 8000 SC/NC + 1 x SKIN	03.00
8002	INCUBATOR 8000 SC/NC + 2 x SKIN	03.00
8003	INCUBATOR 8000 IC	03.00
8004	INCUBATOR 8000 IC + 1 x SKIN	03.00
8005	INCUBATOR 8000 IC + 2 x SKIN	03.00

Available Data

Current Measured Data, Alarm status, Device settings and text messages for Incubator 8000 IC/SC/NC are available from version 1.00.

Commands

Transmitted Commands

Code	Command Specification
30H	Do nothing (NOP)
51H	Initialize Communication (ICC)
52H	Request Device Identification

Processed and responded Commands

Code	Command Specification
24H	Request current DATA
27H	Request current ALARMS
29H	Request current DEVICE SETTINGS
2AH	Request current TEXT MESSAGES
30H	Do nothing (NOP)
51H	Initialize Communication (ICC)
52H	Request Device Identification
55H	Stop Communication

Measured Data

Code	Data Description	Unit	Format
6CH	Air Humidity	%	_XX_
6DH	Air Temperature	°C	XX.X
C3H	Temperature 1	°C	XX.X
BEH	Temperature 2	°C	XX.X
F0H	Inspiratory O2-Concentration	%	XXX_

Temperature 1 = Core temperature

Temperature 2 = Peripheral temperature

Alarm Messages

Air Module

	PRIO 18	Problems with Fan		
CODE	D: LUEFTER INOP	GB: FAN ERR	F: VENT INOP	
CAH	I: ERR VENTOLA	NL: VENT INOP	E: VENT INOP	

	PRIO 20	Ambient Temp. > high Limit		
CODE	D: AMB TEMP "#	GB: AMB TEMP HI	F: AMB TEMP "#	
30H	I: AMB TEMP "#	NL: AMB TEMP "#	E: AMB TEMP "#	

	PRIO 20	Ambient Temperature Sensor inop		
CODE	D: AMB TEMP INOP	GB: AMB TEMP ERR	F: AMB TEMP INOP	
48H	I: AMB TEMP INOP	NL: AMB TEMP INOP	E: AMB TEMP INOP	

	PRIO 8	Ambient Temp. Setting Deviation > 1.5 °C		
CODE	D: AMB TEMP DIF	GB: AMB TEMP DIF	F: AMB TEMP DIF	
6BH	I: AMB TEMP DIF	NL: AMB TEMP DIF	E: AMB TEMP DIF	

	PRIO 3	Lock Setting Temperature > 37 °C active		
CODE	D: AMB TEMP > 37	GB: AMB TEMP > 37	F: AMB TEMP > 37	
4EH	I: AMB TEMP > 37	NL: AMB TEMP > 37	E: AMB TEMP > 37	

Skin Module

	PRIO 20	Temp. 1 – Probe disconnected or fault		
CODE	D: TEMP 1 INOP	GB: TEMP 1 ERR	F: TEMP 1 INOP	
46H	I: ERR TEMP 1	NL: TEMP 1 INOP	E: TEMP 1 INOP	

	PRIO 8	Temp. 1 – Setting Deviation > 0.5 °C		
CODE	D: TEMP 1 DIF	GB: TEMP 1 DIF	F: TEMP 1 DIF	
6DH	I: TEMP 1 DIF	NL: TEMP 1 DIF	E: TEMP 1 DIF	

O2 Module

	PRIO 20	O2 Sensor inoperable		
CODE	D: O2 SENS INOP	GB: O2 SENS ERR	F: CAPT O2 INOP	
43H	I: ERR SENS O2	NL: O2–SENS INOP	E: SENS O2 INOP	

	PRIO 8	O2 Setting Deviation > 5 %		
CODE	D: O2 DIF > 5 %	GB: O2 DIF > 5 %	F: O2 DIF > 5 %	
70H	I: O2 DIF > 5 %	NL: O2 DIF > 5 %	E: O2 DIF > 5 %	

	PRIO 3	Lock O2 Setting > 40 % active		
CODE	D: O2 > 40 %	GB: O2 > 40 %	F: O2 > 40 %	
71H	I: O2 > 40 %	NL: O2 > 40 %	E: O2 > 40 %	

Humidity Module

	PRIO 15	Humidity Sensor inoperable		
CODE	D: FEU SEN INOP	GB: HUM SENS ERR	F: HUM CAP INOP	
34H	I: ERR UMI SENS	NL: HUM SEN INOP	E: HUM SEN INOP	

	PRIO 7	Water Reservoir empty		
CODE	D: WASSERMANGEL	GB: WATER OUT	F: EAU FINI	
9EH	I: ACQUA FINI	NL: WATER STOP	E: AQUA FALTA	

Device Settings

Code	Data Description	Unit	Format
1CH	Air Humidity	%	__XX__
1AH	Air Temperature	°C	_XX.X
1BH	Temperature Skin	°C	_XX.X
01H	Inspiratory O ₂ -Concentration	%	_XXX_

Text Messages

CODE 13H	D: Modul LUFT aktiv	GB: Module AIR active
	F: Module AIR actif	NL: Module LUCHT aktief
	I: Modo ARIA attivo	E: Modo AERE activado


CODE 14H	D: Modul HAUT aktiv	GB: Module SKIN active
	F: Mode PEAU actif	NL: Module HUID aktief
	I: Modo PELLE attivo	E: Modo PIEL activado

CODE 15H	D: Modul O ₂ aktiv	GB: Module O ₂ active
	F: Module O ₂ actif	NL: Module O ₂ aktief
	I: Modo O ₂ attivo	E: Modo O ₂ activado

CODE 16H	D: Modul FEUCHTE aktiv	GB: Module HUM. active
	F: Mode HUMID actif	NL: Module HUMID aktief
	I: Modo HUMID attivo	E: Modo HUMID activado

Specification for Babytherm 8004/8010

Device Connection



Port Specification

Connector	
Type	RS-232-C
Stifte	9-polig, Sub D (female)
	1 Housing
	2 RXD
	3 TXD
	5 GND
Galvanic Isolation	1,5 kV
Position	rear side of Babytherm 8004/8010
	Label: Baby Link®

To connect a PC to Babytherm 8004/8010
the "MEDIBUS-Cable" (Sach-Nr.: 83 06 488) is recommended.

Port-Configuration	
Baudrate	9600 Baud
Databits	8
Startbits	1
Stopbits	1
Parity	none

Gerätekennung

ID Number	Name	MEDIBUS-Version
8006	BABYTHERM 8004	04.00
8007	BABYTHERM 8010	04.00

Available Data

Current Measured Data, Alarm status, Device settings and text messages for Incubator 8000 IC/SC/NC are available from version 1.00.

Commands

Transmitted Commands

Code	Command Specification
30H	Do nothing (NOP)
51H	Initialize Communication (ICC)
52H	Request Device Identification

Processed and responded Commands

Code	Command Specification
24H	Request Current DATA
27H	Request current ALARMS (Codepage 1)
2EH	Request current ALARMS (Codepage 2)
29H	Request current DEVICE SETTINGS
2AH	Request current TEXT MESSAGES
30H	Do nothing (NOP)
51H	Communication initialized (ICC)
52H	Request Device Identification
55H	Stop Communication
4AH	Configure response characteristics
4BH	Set language version

Measured Data

Code	Data Description	Unit	Format
6AH	Mattress temperature	°C	XX.X
6BH	Radiant heater output	%	_XXX
C3H	Temperature 1	°C	XX.X
BEH	Temperature 2	°C	XX.X

Temperature 1 = Core temperature
Temperature 2 = Peripheral temperature

Alarm Messages

Skin Module (Codepage 1)

CODE	PRI0 20	Temp. 1 – Probe disconnected or fault		
	D: TEMP 1 INOP	GB: TEMP 1 ERR	F: TEMP 1 INOP	
	I: ERR TEMP 1	NL: TEMP 1 INOP	E: TEMP 1 INOP	

CODE	PRI0 20	Temp. 2 – Probe disconnected or fault		
	D: TEMP 2 INOP	GB: TEMP 2 ERR	F: TEMP 2 INOP	
	I: ERR TEMP 2	NL: TEMP 2 INOP	E: TEMP 2 INOP	

CODE	PRI0 8	Temp. 1 – Setting Deviation > 0,5 °C		
	D: TEMP 1 DIF	GB: TEMP 1 DIF	F: TEMP 1 DIF	
	I: TEMP 1 DIF	NL: TEMP 1 DIF	E: TEMP 1 DIF	

CODE	PRI0 11	Temp.1 > upper limit		
	D: TEMP 1 "#	GB: TEMP 1 HIGH	F: TEMP 1 "#	
	I: TEMP 1 "#	NL: TEMP 1 "#	E: TEMP 1 "#	

Mattress temperature (Codepage 2)

	PRIO 29	Mattress temperature > 40 °C		
CODE	D: MAT TEMP "#	GB: MAT TEMP HI	F: TEMP MAT "#	
7CH	I: MAT TEMP "#	NL: MAT TEMP "#	E: TEMP MAT "#	

	PRIO 29	Mattress temperature – sensor inoperable		
CODE	D: MT TEMP INOP	GB: MAT TEMP ERR	F: TEMP MT INOP	
7DH	I: ERR TEMP COL	NL: MT TEMP INOP	E: ERR TEMP MAT	

	PRIO 14	Mattress temperature – setting deviation		
CODE	D: MAT TEMP DIF	GB: MAT TEMP DIF	F: TEMP MAT DIF	
7EH	I: TEMP MAT DIF	NL: MAT TEMP DIF	E: TEMP COL DIF	

Radiant heater (Codepage 2)

	PRIO 17	Radiant heater after 15 mins. operation		
CODE	D: RH 15 MIN	GB: RH 15 MIN	F: RH 15 MIN	
7FH	I: RH 15 MIN	NL: RH 15 MIN	E: RH 15 MIN	

Device Settings

Code	Data Description	Unit	Format
1BH	Skin temperature	°C	_XX.X
3EH	Radiant heater output	%	__XXX
3FH	Mattress temperature	°C	_XX.X

Text Messages

CODE 12H	Ventilation mode DS		
	D: Betriebsart DS	GB: Mode DS	F: Mode DS
	NL: Mode DS	I: Modo DS	E: Modo DS


CODE 3BH	Radiant heater is in manual mode		
	D : Strahlungsheizung manueller Mode	GB: Radiant Heater manual Mode	F: Radiateur thermique mode manuel
	NL: Verwarming straler man.-mode	I : Radiatore modo manual	E: Radiador modo manual

CODE 3CH	Radiant heater is in skin temperature control mode		
	D : Strahlungsheizung Haut-Mode	GB: Radiant Heater Skin Mode	F: Radiateur thermique mode peau
	NL: Verwarming straler huid-mode	I : Radiatore modo pelle	E: Radiador modo piel

CODE 3DH	Photo-therapy is on		
	D : Phototherapie an	GB: Photo Therapy on	F: Phototherapie en marche
	NL: Fototherapie aan	I : Fototerapia inserira	E: Fototerapia con

MEDIBUS Specification for Babylog 8000/Babylog 8000 SC

Device Connection



Port Specification

- Connector
- Type RS-232-C
 - 9 pin Sub D (female)
- Pins
- 1 Housing
 - 2 RXD
 - 3 TXD
 - 5 GND
- Galvanic Isolation 1.5 kV
- Location rear side of Babylog 8000 / Babylog 8000 SC
- Label: RS 232
- To connect a PC to Babylog 8000/Babylog 8000 SC the "Medi-Cable" 83 06 488 is recommended.

Port Configuration

- Baudrate 9600 Baud
- Databits 8
- Startbits 1
- Stopbits 1
- Parity none

Device Identification

ID Number	Name	MEDIBUS Version
7000	BABYLOG 8000	03.00
7000	BABYLOG 8000 SC	03.00

Available Data

Data	Software versions Babylog 8000						
	2.00	3.00	3.02	4.02	4.03	4.04	following
Current Measure Data		X	X	X	X	X	X
Device Settings		X	X	X	X	X	X
Realtime Data			X	X	X	X	X
Languages D, GB, F, US, NL		X	X	X	X	X	X
Language E			X			X	X
Languages I, S, J						X	X

Data	Software versions Babylog 8000 SC		
	1.00	1.01	following
Current Measure Data	X	X	X
Device Settings	X	X	X
Realtime Data	X	X	X
Languages D, GB, F, US, NL	X	X	X
Language E		X	X
Languages I, S, J		X	X

The available languages are given for alarm phrases and text messages.

Commands

Transmitted Commands

Code	Command Specification
30H	Do nothing (NOP)
51H	Communication Initialized (ICC)
52H	Request Device Identification

Processed and responded Commands

Code	Command Specification
24H	Request current DATA
27H	Request current ALARMS
29H	Request current DEVICE SETTINGS
2AH	Request current TEXT MESSAGES
30H	Do nothing (NOP)
4AH	Configure Data Response
51H	Initialize Communication (ICC)
52H	Request Device Identification
53H	Request Real time Configuration
55H	Stop Communication

Measured Data

Code	Data Description	Unit	Format
73H	Mean Breathing Press.	mbar	*_XX_
78H	PEEP Breathing Press.	mbar	*_XX_
7DH	Peak Breathing Press.	mbar	*_XX_
80H	Gas Transport Coefficient	mL ² /s	XXXX
85H	Insp. mandatory Tidal Volume	mL	XXX_
86H	High frequency Tidal Volume	mL	XXX_
88H	Tidal Volume	mL	XXX_
B3H	Leakage	%	_XXX
B6H	Spontaneous Fraction Min. Vol.	%	_XXX
B8H	Respiratory Minute Volume (high resolution)	L/min	X.XX
B9H	Respiratory MV (low resolution)	L/min	XX.X
BFH	Insp. mandatory Minute Volume	L/min	X.XX
D6H	Resp. Rate (Vol./Flow)	1/min	XXX_
F0H	Insp. O2	%	XXX_

* Value can be negative. If so, '-' will be the first character.

Realtime Data

Code	Realtime Data	Unit
00H	Airway Pressure	mbar
01H	Flow (insp./exp.)	L/min

Alarm Messages

CODE BEH	PRIOR 8	O2 Measurement inoperable		
	D: FI O2 INOP	GB: % O2 ERR	F: FI O2 INOP	
	I: FI O2 INOP	NL: FI O2 INOP	E: O2 INSP INOP	
	S: % O2 ERR	US: % O2 ERR	J: % O2 ERR	

CODE C1H	PRIOR 8	Flow Measurement inoperable		
	D: FLOW INOP	GB: VOL ERR	F: SPIRO INOP	
	I: FLUSSO INOP	NL: FLOW INOP	E: FLUJO INOP	
	S: VOL ERR	US: VOL ERR	J: VOL ERR	

Device Settings

Code	Data Description	Unit	Format
01H	Insp. Oxygen	%	_XXX_
02H	Inspiratory Flow	L/min	XXX.X
03H	Expiratory Flow	L/min	_XX.X
05H	Inspiratory Time	sec	XX.XX
06H	Expiratory Time	sec	XX.XX
07H	I-Part T _I :T _E	–	_XX.X
08H	E-Part T _I :T _E	–	XXX.X
09H	Frequency IMV (SIMV)	1/min	XXX.X
0BH	PEEP/CPAP	mbar	_XX.X
11H	Apnea Time	sec	_XX.X
13H	Max. insp. Airway Press.	mbar	XXX.X
14H	Trigger Volume	mL	_XX.X
18H	High Minute Volume Limit	L/min	XX.XX
19H	Low Minute Volume Limit	L/min	XX.XX
1DH	Minute Volume Alarm Delay	sec	_XX__
2AH ¹⁾	Frequency HFV	Hz	_XX_
2BH ¹⁾	Amplitude HFV	%	_XXX_
40H ²⁾	set tidal volume VTset	mL	XXX.X
41H ²⁾	Alarm limit panting	bpm	XXX__

¹⁾ Available since device version 4.02 for Babylog 8000.

²⁾ Available since device version 5.00 for Babylog 8000.

Textmessages

CODE	D: Betriebsart IPPV	GB: Mode IPPV	US: Mode CMV	
	F: mode VC	NL: mode IPPV	S: Mode IPPV	
	I: Modo IPPV	E: Modo CMV	J: Mode CMV	
CODE	D: Betriebsart SIPPV	GB: Mode SIPPV	US: Mode A/C	
	F: mode VAC	NL: mode SIPPV	S: Mode SIPPV	
	I: Modo SIPPV	E: Modo A/C	J: Mode A/C	
CODE	D: Betriebsart SIMV	GB: Mode SIMV	US: Mode SIMV	
	F: mode VACI	NL: mode SIMV	S: Mode SIMV	
	I: Modo SIMV	E: Modo SIMV	J: Mode SIMV	
CODE	D: Betriebsart CPAP	GB: Mode CPAP	US: Mode CPAP	
	F: mode VS-PEP	NL: mode CPAP	S: Mode CPAP	
	I: Modo CPAP	E: Modo CPAP	J: Mode CPAP	
CODE	D: Betriebsart DS	GB: Mode DS	US: Mode DS	
	F: mode DS	NL: mode DS	S: Mode DS	
	I: Modo DS	E: Modo DS	J: Mode DS	
CODE	D: Betriebsart VIVE	GB: Mode VIVE	US: Mode VIVE	
	F: mode D.E.V.	NL: mode VIVE	S: Mode VIVE	
	I: Modo VIVE	E: Modo VIVE	J: Mode VIVE	
CODE	D: Betriebsart IPPV + HF	GB: Mode IPPV + HF	US: Mode CMV + HF	1)
	F: mode VC + HF	NL: mode IPPV + HF	S: Mode IPPV + HF	
	I: Modo IPPV + HF	E: Modo CMV + HF	J: Mode CMV + HF	
CODE	D: Betriebsart SIMV+HF	GB: Mode SIMV+HF	US: Mode SIMV+HF	1)
	F: mode VACI+HF	NL: mode SIMV+HF	S: Mode SIMV+HF	
	I: Modo SIMV+HF	E: Modo SIMV+HF	J: Mode SIMV+HF	
CODE	D: Betriebsart CPAP + HF	GB: Mode CPAP + HF	US: Mode CPAP + HF	1)
	F: mode VS-PEP + HF	NL: mode CPAP + HF	S: Mode CPAP + HF	
	I: Modo CPAP + HF	E: Modo CPAP + HF	J: Mode CPAP + HF	
CODE	D: Betriebsart PSV	GB: Mode PSV	US: Mode PSV	
	F: mode AI	NL: mode PSV	S: Mode PSV	
	I: Modo PSV	E: Modo PSV	J: Mode PSV	
CODE	D: Betriebsart PSV+VG	GB: Mode PSV+VG	US: Mode PSV+VG	
	F: mode AI+VG	NL: mode PSV+VG	S: Mode PSV+VG	
	I: Modo PSV+VG	E: Modo PSV+VG	J: Mode PSV+VG	
CODE	D: Betriebsart SIMV+VG	GB: Mode SIMV+VG	US: Mode SIMV+VG	
	F: mode VACI+VG	NL: mode SIMV+VG	S: Mode SIMV+VG	
	I: Modo SIMV+VG	E: Modo SIMV+VG	J: Mode SIMV+VG	
CODE	D: Betriebsart SIPPV+VG	GB: Mode SIPPV+VG	US: Mode A/C+VG	
	F: mode VAC+VG	NL: mode SIPPV+VG	S: Mode SIPPV+VG	
	I: Modo SIPPV+VG	E: Modo A/C+VG	J: Mode A/C+VG	

1) Available for Babylog 8000 since device version 4.02



Directive 93/42/EEC
concerning Medical Devices

Dräger Medizintechnik GmbH

Germany

🏠 Moislinger Allee 53 – 55
D-23542 Lübeck

☎ +49 451 8 82 - 0

📠 26 80 70

FAX +49 451 8 82-20 80

💻 <http://www.draeger.com>